



TIN HANG TECHNOLOGY LIMITED

Units 1904-6, 19/F., No. 340 Queen's Road Central, Hong Kong

Tel: (852)-2817 2121 Fax: (852)-2580 7763

Website: <http://www.tinhangtech.com> Email: sales@tinhangtech.com

PROMOTION

Whatman® Leadership in separations technology for the life sciences

We are proud to announce that Tin Hang Technology Limited is now the

Exclusive Distributor of **Whatman Schleicher & Schuell** in Hong Kong.

Whatman S&S is a global leader in separations technology and is known throughout the scientific community for providing innovative products and solutions. She offers comprehensive products in several areas: Filter Papers and Membrane, Filtration Devices, Multiwell Plates, Nucleic acid Sample Preparation & Analysis, Extraction and Chromatography.

ALL Filter Papers, pH Indicator/paper and Benchkote Sheets

Buy 10 get 2 for free!



13mm, 25mm Syringe Filters LESS 10% Discount

(Price list on the last page)

There are free samples available. If you are interested, please do not hesitate to contact us!


~Promotional Period until End of June~

Blotting Membranes:

- PROTRAN®**
- Universal blotting membrane
 - 100% pure nitrocellulose, no cellulose acetate
 - High binding, low background

- NYRRAN®**
- Positively charged nylon membrane
 - Better binding of negatively charged DNA/RNA

- WESTRAN®**
- Clear Signal PVDF membrane
 - Protein blotting membrane



Blotting Paper:

- 3MM CHR**
- Pure cellulose produced from highest quality cotton linters
 - Manufactured and tested for chromatographic techniques



~There are different dimensions for your needs!~

Blotting Devices:

- TurboBlotter™**
- Rapid (**1 hr**) downward DNA, RNA Blotting system
 - Ready-to-blot kit
 - Easy to set-up, NO heavy weights or vacuum required
 - Economical design uses **less buffer (~125-200ml)** and blotting paper



All the above Blotting Products


Less 10% Discount!

~Promotional Period until End of June~

Nucleic Acid Sample Preparation and Analysis:

FTA[®] Card, FTA[®] Elute

- For the Collection, Transport, Storage and Rapid Purification of Nucleic Acids
- Simple, Rapid nucleic acids capture
- Store for years at room temperature
- Fast way to isolate DNA for PCR!



Related Products:

FTA Reagent

- Removes heme, PCR inhibitors, PCR inhibitors and other contaminants
- Ensures superior quality DNA for analysis

FTA Kit

- 25 FTA Micro Card, 2 Purification reagent (25ml), 2 Harris Micro Punches with cutting mat and instruction.

~Less 10% Promotional Discount until End of June~

Interested? Please do not hesitate to contact us!



天恒科技有限公司

TIN HANG TECHNOLOGY LIMITED

香港皇后大道中340號19樓1904-1906室

Units 1904-1906, 19/F, No. 340 Queen's Road Central, Hong Kong

Tel: (852) 2817 2121 Fax: (852) 2580 7763

E-mail: sales@tinhangtech.com

Web site: www.tinhangtech.com