

# Steritop/Stericup Filter Units

## Ideally suited for sterile filtration of cell culture media, buffers & reagents



### Sterile & ready-to-use

### Faster Filtration

Millipore Express Plus membrane can achieve up to 50% faster than competitive 0.2 µm PES (Figure 1)

### Low Binding

Membrane's proprietary hydrophilic surface ensures that **less protein binding** (Figure 2)




Figure 1 47 mm membrane disks with 0.22 µm pore size were challenged with DMEM with 10% FBS (500 mL)




Figure 2 Membrane disks with a 0.22 µm pore size were offered a 1 mg/mL solution of <sup>125</sup>I labeled IgG. The chart shows protein binding after incubation

### Please contact:

#### Millipore China Limited

Room 1701, 17/F., Harcourt House, 39 Gloucester Road, Wanchai, Hong Kong

Tel : (852) 2803 9265 Fax : (852) 2513 0313 Web site : <http://www.millipore.com/sterile>

or



Millipore Distributor:

#### Tin Hang Technology Limited

Unit 1904-1906, 19/F., No.340 Queen's Road Central, Hong Kong

Tel: (852) 2817 2121 Fax: (852) 2580 7763 E-mail: [sales@tinhangtech.com](mailto:sales@tinhangtech.com)  
Web site : <http://www.tinhangtech.com>



## Stericap™ Filter Units Vacuum-driven filtration and storage devices

Units with fast flow and low protein binding 0.22 µm Millipore (PES) membrane

for sterilization of aqueous solution; supplied sterile.



Product Name	Funnel Capacity, mL	Receiver	Qty/pk	Catalogue No.
Stericap-GP	150	150mL bottle with cap	12	SCGPU01RE
	250	250mL bottle with cap	12	SCGPU02RE
	500	500mL bottle with cap	12	SCGPU05RE
	500	1,000mL bottle with cap	12	SCGPU10RE
	1,000	1,000mL bottle with cap	12	SCGPU11RE



## Steritop™ Filter Units Vacuum-driven bottle-top devices

Units with fast flow and low protein binding 0.22 µm Millipore (PES) membrane

for sterilization of aqueous solution; supplied sterile.



Product Name	Funnel Capacity, mL	Receiver	Qty/pk	Catalogue No.
Steritop-GP	150	33 mm thread	12	SCGPS01RE
	250	33 mm thread	12	SCGPS02RE
	500	33 mm thread	12	SCGPS05RE
Steritop-GP	150	45 mm thread	12	SCGPT01RE
	250	45 mm thread	12	SCGPT02RE
	500	45 mm thread	12	SCGPT05RE

# “Buy more Get more Discount”

**3 Packs Less 10%**

**5 Packs Less 20%**

*Promotion Until December, 2008*

### Order Form

Date:	_____
Name:	_____
University/Hospital:	_____
Department:	_____
Location:	_____
Phone:	_____
Fax:	_____
E-mail:	_____

Qty/pk	Cat. No.	Qty/pk	Cat. No.
_____	X SCGPU01RE	_____	X SCGPS01RE
_____	X SCGPU02RE	_____	X SCGPS02RE
_____	X SCGPU05RE	_____	X SCGPS05RE
_____	X SCGPU10RE	_____	X SCGPT01RE
_____	X SCGPU11RE	_____	X SCGPT02RE
_____		_____	X SCGPT05RE

Please send The Following

Fax: 2580 7763